

GlideScope AVL Reusable

VIDEO LARYNGOSCOPE

When time is critical, make your go-to device the GlideScope® AVL.

Ideal for demanding OR, ED, and ICU settings, the AVL provides a consistently clear view of the airway, and enhanced maneuverability enabling fast intubations.

The AVL system gives you a choice of 4 reusable blades to cover a wide range of patients—from infant to morbidly obese.

Add real-time recording and an onboard tutorial, and you'll see why GlideScope AVL is the one you go to for advanced video laryngoscopy.

AVL=ADVANCED AIRWAY VIEWS


- Digital color monitor with DVD-quality resolution
- ► Reusable video laryngoscope with digital-resolution camera
- ► Innovative blade sizing, angulation, and camera positioning
- ► Reveal[™] anti-fog feature, with a rapid heating profile to resist lens clouding and secretions

AVL=ADVANCED USER FEATURES

- Integrated, real-time recording helps confirm tube placement and facilitate teaching
- ► Onboard 4-Step Technique tutorial makes AVL easy to use, learn, and teach
- ► Ergonomic design and rugged construction for multiple clinical applications
- ► Reusable blades in 4 sizes: GVL® 2, 3, 4, 5
- ► Also available in 6 single-use sizes



Advanced Video Laryngoscopy: Real-time recording. DVD clarity.





GlideScope AVL Reusable

VIDEO LARYNGOSCOPE

Applications

- First-use intubations, replacing direct laryngoscopy
- Normal or restricted oropharyngeal views
- · Cormack-Lehane Grades I-IV laryngeal views
- Trauma airways—excellent when dealing with blood and secretions in the airway
- Nasal tracheal intubation
- Video-guided tube exchange in the ICU


System Components

- Digital, DVD-quality color monitor
- Reusable video laryngoscope with digital camera
- Reusable AVL blade in a choice of four (GVL®) sizes
- GlideRite® Rigid Stylet
- Choice of mobile stand or IV pole mount
- In-Service CD including User's Manual and Quick Reference Cards

GlideScope® Products and Systems: Designed for 1st Pass Success

REUSABLE SYSTEM			
GVL 2	GVL 3	GVL 4	GVL 5
4 - 20 kg	10 kg - Adult	40 kg - Morbidly Obese	40 kg - Morbidly Obese
Length (blade tip to handle): 51 mm Height at camera: 11 mm Width at camera: 14 mm	Length (blade tip to handle): 82 mm Height at camera: 16 mm Width at camera: 20 mm	Length (blade tip to handle): 103 mm Height at camera: 14 mm Width at camera: 27 mm	Length (blade tip to handle): 103 mm Height at camera: 14 mm Width at camera: 27 mm





Note: System configurations vary.

*The GVL 5 is designed to accommodate anatomical anomalies sometimes associated with bariatric patients.

Υ...

Verathon Medical Europe B.V. Linnaeusweg 11 / 3401 MS IJsselstein / The Netherlands Tel. +31 (0)30 68 70 570 / Fax +31 (0)30 68 70 512

GlideScope, the GlideScope symbol, GVL, GlideRite, Verathon, and the Verathon Torch symbol are trademarks of Verathon Inc. © 2011, 2012 Verathon Inc.

GlideScope® video laryngoscope systems are CE marked in accordance with the Medical Device Directive, and the Verathon Inc. quality system is Quality System Certified to ISO 13485:2003 standards.